

Vdrive

The ultimate in performance veterinary power tools



deSoutter
MEDICAL

V DRIVE - MODULAR SYSTEM



Order code	Model	Description	Speed rpm
1291774	V-MBQ-707	Single Trigger	0 - 1350
1291884	V-MBQ-708	Twin Trigger	0 - 1350

ATTACHMENTS

	Order code	Model	Description	Speed rpm	Cannulation Ø mm
	18320	V-DQ-708	0.5 - 6.4 mm (1/4") keyed	0 - 1350	4.4
	18330	V-DQ-708	Small A.O. Synthes/ASIF	0 - 1350	2.1
	18440	V-RQ-708	0.5 - 6.4 mm (1/4") keyed	0 - 330	4.4
	18290	V-KQ-707	Sagittal saw Uses S89 type blades	0 - 13500 cpm	
	18300	V-NQ-707	Small sagittal saw Uses S88 type blades	0 - 20000 cpm	
	18270	V-OQ-708	TPLO oscillating saw Uses Slocum blades	0 - 13500 cpm	
	18690	V-OQ-708	TPLO oscillating saw Synthes	0 - 13500 cpm	
	18310	V-WQ-707	0.7 - 4.0 mm wire capacity	0 - 1350 cpm	4.0

STERILISATION CONTAINERS



STANDARD CONTAINERS AND UNIVERSAL TRAY

Universal sterilisation case and tray order code: 9350
Base order code: 630833
Lid order code: 630843



Wire sterilisation container and lid.
Order code: 655073

Universal tray order code: 631343
Blue mat order code: 616153

BATTERY AND ELECTRICAL COMPONENTS



The BC-700 Battery Charger requires minimal user input, simply insert battery and the BC-700 will charge, recover and recondition it to its optimum state.

BC-700 MULTI BAY CHARGER AND MODULE OPTIONS

Model	Charger module description	Charger module - Order code	For battery type
15700	Without charger module fitted	-	-
16210	With 4-off CA-701 charger modules	15760 	AB-701/AB-702
16220	With 4-off CA-700 charger modules	15750 	SB-703/SB-704
16230	With 4-off CA-703 charger modules	18120 	AB-703/AB-704

Mains power cordsets -
3 metres long.
Please specify - Order code

UK	8290
Euro	8300
US/Canada	8310
Australia/New Zealand	8320
Swiss	14700

The BC-700 charger modules can be quickly and easily changed to suit the batteries being charged.



SINGLE BAY CHARGERS

BC-705	For battery type AB-701/AB-702	Order code: 17980
BC-706	For battery type SB-703/SB-704	Order code: 17990
BC-707	For battery type AB-703	Order code: 18000
BC-708	For battery type AB-704	Order code: 18010

LITHIUM ION BATTERIES

Why choose Lithium Ion battery technology?

De Soutter Medical Lithium Ion batteries offer the surgical team and sterile services a number of tangible benefits over conventional battery technology, including:

- Li-Ion technology improves battery 'capacity to weight' ratio by around 100% when compared to NiMH or NiCad chemistries.
- Significantly extends the useable run time of the power tool.
- Eliminates 'memory effect' normally associated with other battery types.
- De Soutter Li-Ion technology readily withstands repeated machine washing and autoclaving.
- Unique **fuel gauge** incorporated within De Soutter Li-Ion batteries allows the surgical team to instantly assess the battery capacity - both before and during use.

Sterile batteries - Li Ion				
Model	Description	Order code	Voltage	Capacity mAh
SB-703	Small battery	17210	13.2V	1100
SB-704	Large battery	17220	13.2V	2500

Aseptic batteries					Battery housing	Order code	Battery shield	Order code
AB-701	Small battery - NiMH	15340	10.8V	1950	V-AH-701	18340	AS-701	15890
AB-702	Large battery - NiMH	15350	14.4V	1950	V-AH-702	18350	AS-702	15900
AB-703	Small battery - Li Ion	17910	10.95V	2000	V-AH-703	18080	AS-703	18100



Vdrive

UNITED KINGDOM (INTERNATIONAL HEAD OFFICE)

De Soutter Medical Limited
Halton Brook Business Park, Weston Road
Aston Clinton, Aylesbury, Bucks HP22 5WF
Tel: +44 (0) 1296 634 000
Fax: +44 (0) 1296 634 033
Email: info@de-soutter.com

GERMANY

De Soutter Medical Germany
Niederlassung Deutschland, Kladenfloss,
D-66625, Nohfelden
Tel: +49 (0) 68 52-99 12 46
Fax: +49 (0) 68 52-99 12 47
Email: deutschland@de-soutter.com

AUSTRALIA

De Soutter Medical Australia Pty Ltd
2/12-14 Apollo Drive, Hallam
Victoria 3803
Tel: +61 (0) 3 9702 4441
Fax: +61 (0) 3 9702 4484
Email: australia@de-soutter.com

ITALY

De Soutter Medical Italia
Località Fornace SNC
27022 Casorate Primo - PV
Tel: +39 (0) 2 9009 4098
Fax: +39 (0) 2 9009 2673
Email: italy@de-soutter.com

AUSTRIA

De Soutter Medical Austria
Zweigniederlassung Österreich
Dietrichsteingasse 10, A-3400 Klosterneuburg
Tel: +43 (0) 676 96 71 770
Fax: +43 (0) 2243 21 656
Email: austria@de-soutter.com

NETHERLANDS

De Soutter Medical Netherlands
Gelderlandhaven 2X
3433 PG Nieuwegein
Nederland
Tel: +31 (0) 850491480
Fax: +31 (0) 850491489
Email: netherlands@de-soutter.com

BELGIUM

De Soutter Medical Belgium
In De Bruel 30
3620 Lanaken
Tel: +32 (0) 89/47 15 37
Fax: +32 (0) 89/70 12 19
Email: belgium@de-soutter.com

UNITED STATES OF AMERICA

De Soutter Medical USA, Inc.
224 Rolling Hill Road, Suite 12A
Mooresville, NC 28117
Tel: +1 (704) 655 9040
Fax: +1 (704) 987 2035
Email: usa@de-soutter.com

FRANCE

De Soutter Medical France
949 Avenue du Parc des Expositions
33260 La Teste de Buch
Tel: +33 (0) 5 56 54 89 36
Fax: +33 (0) 9 70 61 37 60
Email: france@de-soutter.com

